

Goodman HTS Pipeline User Manual

version 0.1

Simon Torres

July 04, 2017

Contents

| | |
|--|----------|
| Goodman High Throughput Spectrograph Pipeline - User Manual | 1 |
| Testing the Pipeline | 1 |
| Short Pipeline Overview | 1 |
| VPN Connection | 1 |
| VNC Connection | 1 |
| Using the Terminal | 1 |
| Using Graphical VNC Clients | 1 |
| Getting Help | 2 |
| Running the Pipeline | 3 |
| redccd | 3 |
| redspec | 3 |
| Overview | 3 |
| Goodman HTS | 3 |
| What to expect from this Pipeline | 4 |
| Getting the Data | 4 |
| Package Content | 4 |
| Installing | 4 |
| Install prerequisites | 4 |
| Get Goodman HTS Pipeline | 4 |
| Using built-in Help | 4 |
| Importing Parts of the Code | 4 |

Goodman High Throughput Spectrograph Pipeline - User Manual

This document is the User Manual for Goodman HTS Pipeline. It will provide a quick overview of the main features of the pipeline, installation instructions, restrictions and guidelines for using it.

Testing the Pipeline

We have set up a server that you can access using VNC. You can test it on your own data or the data we have provided for test. Also, you can test it on your own machine but we can't give you support, although we provide full documentation for installing and testing it.

Short Pipeline Overview

The Goodman Pipeline is distributed as a single package but the full process is split in two: `redccd` and `redspec`. The first does the 2D image reduction, going from raw data to flatfielded data, and the second (`redspec`) takes a 2D image containing one or more spectrum and creates a wavelength calibrated 1D spectrum.

VPN Connection

If you are connecting from outside AURA you must connect through the VPN if you don't have the information, ask your *support scientist* to provide it for you.

VNC Connection

The server is `soardata3` IP Address `139.229.15.174` and the vnc password is `vnc4goodmanpipe`.

Using the Terminal


Open a terminal, and assuming you have installed `vncviewer`.

```
vncviewer soardata3:1
```

You will be asked to type in the *password* provided above.

Using Graphical VNC Clients

Using a graphical VNC client is quite similar and intuitive



In this case the *IP address* was used, which is equivalent and sometimes better.

Getting Help

The Goodman Pipeline comes with full documentation and you can get help by using the `--help` argument. For instance:

For `redccd`:

```
bash$ redccd --help
usage: redccd [-h] [--cosmic] [--ignore-bias] [--auto-clean]
             [--saturation <Value>] [--raw-path raw_path]
             [--red-path red_path] [--debug] [--log-to-file]
             [--flat-normalize <Normalization Method>]
             [--flat-norm-order <Order>] [--dcr-par-dir <dcr.par directory>]
             [--keep-cosmic-files]
```

Goodman CCD Reduction - CCDreductions for Goodman spectroscopic data

optional arguments:

| | |
|--|--|
| <code>-h, --help</code> | show this help message and exit |
| <code>--cosmic</code> | Clean cosmic rays from science data. |
| <code>--ignore-bias</code> | Ignore bias correction |
| <code>--auto-clean</code> | Automatically clean reduced data directory |
| <code>--saturation <Value></code> | Saturation limit. Default to 55.000 ADU (counts) |
| <code>--raw-path raw_path</code> | Path to raw data. |
| <code>--red-path red_path</code> | Path to reduced data. |
| <code>--debug</code> | Show detailed information of the process. |
| <code>--log-to-file</code> | Write log to a file. |
| <code>--flat-normalize <Normalization Method></code> | Choose a method to normalize the master flat for spectroscopy. Choices are: mean, simple (model) and full (fits model to each line). |
| <code>--flat-norm-order <Order></code> | Defines the order of the model to be fitted. |
| <code>--dcr-par-dir <dcr.par directory></code> | Directory of default dcr.par file. |
| <code>--keep-cosmic-files</code> | After cleaning cosmic rays with dcr, do not remove the input file and the cosmic rays file. |

And for `redspec`:

```
bash$ redspec --help
usage: redspec [-h] [--data-path <Source Path>]
              [--proc-path <Destination Path>]
              [--search-pattern <Search Pattern>]
              [--output-prefix <Out Prefix>] [--extraction <Extraction Type>]
              [--reference-files <Reference Dir>] [--interactive] [--debug]
              [--log-to-file] [--save-plots] [--plot-results]
```

Extracts goodman spectra and does wavelength calibration.

optional arguments:

| | |
|--|---|
| <code>-h, --help</code> | show this help message and exit |
| <code>--data-path <Source Path></code> | Path for location of raw data. Default <./> |
| <code>--proc-path <Destination Path></code> | Path for destination of processed data. Default <./> |
| <code>--search-pattern <Search Pattern></code> | Pattern for matching the goodman's reduced data. |
| <code>--output-prefix <Out Prefix></code> | Prefix to add to calibrated spectrum. |
| <code>--extraction <Extraction Type></code> | Choose a which extraction to perform. Simple is a sum |

```

                                across the spatial direction after the background has
                                been removed. Optimal is a more advanced method that
                                considers weights and profilefitting.
--reference-files <Reference Dir>
                                Directory of Reference files location
--interactive
                                Interactive wavelength solution.Disabled by default.
--debug
                                Debugging Mode
--log-to-file
                                Write log to a file
--save-plots
                                Save all plots in a directory
--plot-results
                                Show wavelength calibrated spectrum at the end.

```

Running the Pipeline

The pipeline is designed to work even if you put no arguments but this not always the best. For well behaved data this might be useful or if you have a large survey with data obtained in an ordered and systematic way.

redccd

For redccd I suggest using `--cosmic` and `auto-clean` also you might want to consider `--saturation <new value>` to change the saturation level if you get all your flats rejected due to saturation. Sometimes there is a hot column at the end that produced very high values.

```
redccd --cosmic --auto-clean
```

In case you want to use `--saturation` here is an example:

```
redccd --cosmic --auto-clean --saturation 70000
```

This changes the saturation level to *70000 ADU* in this context the saturation value works as a threshold for rejecting images.

By default, redccd puts reduced data in a subdirectory `RED`, you can provide a different one by using `--red-path`.

An image `image_file.fits` that has been fully (and properly) processed should have the new name (including the reduced data folder):

```
cfzsto_image_file.fits
```

Where *c* stands for *cosmic ray rejected*, *f* for flatfielded, *z* for zero or bias corrected, *s* for slit trimmed, *t* for trimmed and *o* for overscan corrected.

redspec

By default redspec will search for images with the prefix `cfzsto` in case you have produced a different prefix you can change it by using `--search-pattern`

You can just run redspec in case everything is the default but if this is the first time you run the pipeline I suggest:

```
redspec --plot-results
```

In that way two important plots will be shown full screen, the comparison lamp fitted to a reference comparison lamp and some values for the wavelength solution fit and the extracted spectrum plotted with the wavelength solution.

Overview

Goodman HTS

The Goodman High Throughput Spectrograph was developed by the [Goodman Laboratory](#) at University of North Carolina and is currently in operation at the [SOAR Telescope](#) is an *imaging spectrograph* which means you can obtain Images and Spectra. It is a highly customizable Instrument, has two cameras, *Red* and *Blue*. As you might have guess already the *Red* Camera

What to expect from this Pipeline

Getting the Data

Package Content

Installing

Install prerequisites

Get Goodman HTS Pipeline

Using built-in Help

Importing Parts of the Code

You can import almost every function in this pipeline:

```
1 from goodman_ccd.core import spectroscopic_extraction
2 # do some stuff here
3 print('Hello World')
```